



Master Complete

INTRODUCTION

- Laser engraved box
- Stainless steel case
- Embedded thermal printer
- Water and dust resistant keyboard
- User-friendly interface
- Serial port Rs232
- Analysis of 3 types of milk
- No thermal or mechanical treatment of the sample before analysis
- Low power consumption
- Speed - 60sec.
- Autocalibration (No PC required)
- No use of chemical reagents
- Low Cost



TECHNICAL DATA	
Power supply	AC 95V - 250V 50-60Hz
Input for use in car	DC 12V - 14V
Product dimension	235/235/115
Weight	3 kg.

PARAMETERS

- Fatness (FAT)
- Solids non fat (SNF)
- Milk density
- Protein
- Added water to milk
- Freezing point
- Lactose
- Solids

We offer to your attention the offer of Milkotester - Ltd for the design, manufacture and service of ultrasonic milk analyzers and laboratory equipment. Milk analyzers Master Series are created using multiple developments and solutions for achieving the high accuracy of measurement indicators. Products manufactured by Milkotester have proven their quality in over 50 countries.

Bulgaria-4470 Belovo
Hristo Botev No49 Str.
www.milkotester.com